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ENGLISH

Attention

Please read carefully the information contained in this folder.

Inexpert handling and care, as well as improper use may result in premature wear or hazards for patients and users.

Use for the intended purpose

Cables for electrosurgery are designed either to conduct electrical power from the output of a high-frequency generator to the instrument or to connect a neutral electrode with the generator.

Attention: Instruments and cables for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.

Use and safety instructions

- All cables have to be completely cleaned, disinfected and sterilised before initial use and any other use.
- The cables have to be submitted to a visual inspection and a functional test before each use.
- Never use damaged cables.
- Store the cables in a clean, cool and dry place. Protect them from damage.
- Handle the cables with care, and do not pinch them.
- Make sure that the cable is correctly inserted in the socket or connected with the plug.

Reprocessing

Preparation and transport

Remove coarse dirt from the cables immediately after each use. Do not use fixation agents or hot water (>40 °C).

Machine reprocessing

Cleaning

Place the cables in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. Prerinse for 1 min. with cold water
2. Discharging
3. Prerinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55 °C with a 0.5% alkaline or at 45 °C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40 °C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40 °C).
10. Discharging

Disinfection

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying

Dry the cables by carrying out a drying cycle of the cleaning/disinfection machine.

If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing

Cleaning

Prepare a cleaning bath according to the manufacturer's instructions.

1. Rinse products with cold tap water (<40 °C): until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submersed. Observe residence time according to the manufacturer's instructions.
3. Clean the cable in the bath manually using a soft brush. Brush all surfaces several times.
4. Rinse the products thoroughly with running tap water to remove the cleaning agents without residue.

Disinfection

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the cables in the disinfectant bath and observe the specified residence time. Rinse the cables thoroughly with fully demineralised water to remove the disinfectant without residue.

Drying

Dry manually by using a lintfree cloth and, in particular, for drying cavities, sterile compressed air.

Functional test and packaging

Perform visual inspection for cleanliness. If necessary, repeat the reprocessing process until the instrument is optically clean. Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.

- 3 pre-vacuum phases with a pressure of at least 60 mbar.

- Heating up to a sterilisation temperature of min. 132 °C and max. 137 °C
- Shortest exposure time: 3 min.
- Drying time: at least 10 min.

Validation of the reprocessing

Information on the validation of the reprocessing and detailed handling instructions can be downloaded from www.bissinger.com.

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