Distributed by:

Black & Black Surgical, Inc.
5175 South Royal Atlanta Dr.
Tucker, GA 30084 USA
Tel: 877-252-2517 (in USA)
770-414-4880 (outside USA)
Fax: 770-414-4879

HANDHABUNGSHINWEIS

INSTRUCTIONS FOR USE

INSTRUCTIONS DE SERVICE

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Bipolare Pinzetten

Bipolar Forceps

Pinces bipolaires

Pinze bipolari

Pinzas bipolares

Black & Black Surgical
B 68210-06 - B 68711-1

CE 0297

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Reinigungsmittel (maschinell):

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Reinigungsmittel (manuell):

Enzol Enzym. Detergent, Johnson&Johnson

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Bipolar forceps

Attention

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

Intended Use

Bipolar forceps are used for grasping, dissection and coagulation of tissues. They have to be connected, by means of an appropriate bipolar cable, to the bipolar output of a high-frequency generator and may only be used with bipolar coagulation current.

Maximum output voltage of the generator U_{max} :

500 V_r or 300 V_r (depending on the type, see catalogue)

Appropriate connecting cables:

Bissinger bipolar cable / flat plug, item no. 801 00xxx.

Bissinger bipolar cable / 2-pin plug, item no. 801 00xxx.

Attention: Instruments for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.

Contraindications

Incidents which have been reported in connection with the use of bipolar systems:

Unintended activation with resulting tissue injury on the wrong spot and/ or damage to the equipment.

Fire in connection with surgical drapes and other inflammable materials.

Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.

Explosions caused by sparks in the proximity of inflammable gases.

Perforation of organs. Sudden severe bleedings.

Bipolar forceps have proved inefficient for tubular sterilisation or coagulation in the context of sterilisation and should therefore not be used for this purpose.

Use and safety instructions

The non-observance of the present use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.

- Before initial use and any other use, all instruments have to be completely cleaned, disinfected and sterilised and their functioning has to be checked.

- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas, such as blades, tips, notches, locking and blocking devices as well as all mobile parts, insulations and ceramic elements have to be checked carefully.

- Never use damaged instruments.

- Never use the instruments in the presence of flammable or explosive substances.

- The instrument may not be laid down on the patient.

- Coagulation should only be performed if the contact surfaces are visible and ensure a good contact to the tissue selected for coagulation. Do not touch any other metallic instruments, trocar sleeves, optics or things like that during use.

Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

Reprocessing

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling.

Instrument for electrosurgery are by nature subject to increased wear depending on the type and time of use.

Preparation and transport

Remove coarse dirt from the instruments immediately after each use. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container.

Machine reprocessing

Cleaning

Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.

1. Pre-rinse for 1 min. with cold water
2. Discharging
3. Pre-rinse for 3 min. with cold water
4. Discharging
5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
6. Discharging
7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
8. Discharging
9. Rinse for 2 min. with warm tap water (>40°C).

10. Discharging

Disinfection

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.

If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing

Ultrasonic Pre-Cleaning

1. The instrument must be inserted in an ultrasonic bath with 0.5% enzymatic cleaning detergent. Ultrasonic must be applied for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse completely with cold water to remove the cleaning detergent.

Cleaning

Prepare a cleaning bath according to the manufacturer's instructions.

1. Rinse products with cold tap water (<40°C) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
2. Place products in the prepared cleaning bath so that they are completely submerged. Observe residence time according to the manufacturer's instructions.
3. Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
4. **The following step only applies to channels and the insides of tubes:** Push the brush into and out of the tubes at least six times. Rinse the tubes with DI water. Repeat the procedure.
5. Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

Disinfection

Prepare a disinfectant bath according to the instructions of the disinfectant manufacturer. Place the instruments in the disinfectant bath and observe the specified residence time. Rinse the products thoroughly with fully demineralised water to remove the disinfectant without residue.

Drying

Manual drying is carried out using a lintfree cloth and, in particular, for drying cavities and channels, sterile compressed air.

Functional test and packaging

Perform visual inspection for cleanliness; if required, perform an assembly and functional test according to the operating instructions.

If necessary, repeat the reprocessing process until the instrument is optically clean. Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

CLARIS non-stick bipolar forceps

The polished noble metal forceps tips of the CLARIS bipolar forceps may tarnish like silver. This does not lead to any malfunction

Bipolar forceps with irrigation

The wire insert comprised in the scope of supply should be always inserted in the irrigation channel - except during utilisation and cleaning - in order to prevent any clogging. The irrigation channel has to be rinsed very thoroughly during cleaning. The passage has to be checked after cleaning.

Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.

- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C
- Shortest exposure time: 3 min.
- Drying time: at least 10 min.

Storage

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

Repairs

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. For any question on these

matters, please contact either the manufacturer or your medico-technical department.

Attention: Defect products must pass the complete reprocessing process before being returned for repair.

Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):

Neodisher FA by Dr. Weigert (alkaline)

Endozime by Ruhof (enzymatic)

Cleaning agents (manual cleaning):

Enzol Enzym, detergent by Johnson&Johnson

Disinfectants (manual disinfection):

Cidex OPA, Johnson&Johnson

Neutrlising agent:

Neodisher Z by Dr. Weigert

Cleaning and disinfection device:

Miele G 7736 CD

Miele insert module E 327-06

Miele MIS module E 450

For details, see report.

SMP GmbH # 01707011901 (machine cleaning)

MDS GmbH # 135196-10 (man. cleaning/disinfection)

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If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care.

This applies particularly to blades, fine tips and other sensitive areas.

Warranty

Günter Bissinger Medizintechnik GmbH exclusively supplies tested and faultless products to the customers. All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexpert way.

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SMP GmbH # 01707011901 (limpieza mecánica)
MDS GmbH # 135196-10 (limpieza/desinfección manual)
MDS GmbH Testbericht 084183-10 (esterilización)
DQS Medizinprodukte GmbH