**Distributed by:**  
**Black & Black Surgical, Inc.**  
 5175 South Royal Atlanta Dr.  
 Tucker, Ga 30345. USA  
 Tel: 877-252-2517 (in USA)  
 770-414-4880(outside USA)  
 Fax: 770-414-4879

**HANDHABUNGSHINWEISE**  
**HANDLING INSTRUCTIONS**  
**INSTRUCTIONS DE SERVICE**  
**ISTRUZIONI PER L'USO**

**BiTech Bipolare Scheren**  
**BiTech Bipolar Scissors**  
**Ciseaux bipolaires BiTech**  
**Forbici bipolari BiTech**

Artikel-Nr. B68750 - B68751  
 Article No. B68780 - B68781  
 N° d'article  
 No. di articolo

Bissinger Medizintechnik GmbH  
 Hans-Theisen-Str.1  
 79331 Teningen / Germany

Tel.: ++49 7641 9 14 33 0  
 Fax: ++49 7641 9 14 33 33  
 Email: info@bissinger.com

**Deutsch**  
**BiTech Bipolare Scheren**

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**English**  
**BiTech Bipolar Scissors**

**Intended Use**  
 BiTech bipolar scissors are designed for dissecting, cutting and bipolar coagulation of tissue. They have to be connected by a suitable bipolar cable to the appropriate output of a high frequency electrosurgical generator. The coagulation current is activated by means of a footswitch. The max. output voltage of the generator may not exceed 300 Vp. The output electrical capacity should be between 40 and 70 watt (coagulation mode).  
 Suitable bipolar cables: Bissinger Art.No. 801 00081 – 801 00088.

**Note:** Instruments for electrosurgery should only be used by individuals who have been specially trained in their use.

**Safety Instructions**  
 Examine the instrument for damage before each use. Bipolar scissors have high-quality ceramic parts which have to be treated with particular care and protected against breakage.  
 Do not use in the proximity of flammable or explosive materials.  
 Do not use on patients with pacemakers.  
 Only activate power when the contact area is in full view.  
 Do not come into contact with other metal instruments.  
 Never attempt to make repairs yourself.

**Initial use**  
 All new instruments must be cleaned and sterilized before initial use.  
 Carry out a visual inspection and test functions before each use. Clean and sterilize before each use.

**Storage and Transportation**  
 Instruments must be stored in a clean, cool and dry area. Protect from mechanical damage.  
 Handle with extreme care, do not throw or drop !

**Cleaning and Disinfection**  
 Instruments should be thoroughly cleaned immediately after each use to prevent the adhesion of particles and dried secretions.  
 Use commercially available instrument cleaning and disinfection agents. Agents with neutral pH (7.0) are best. It is essential to follow the manufacturer's instructions. Instruments can be cleaned and disinfected in the instrument washing machine up to +95°C (203°F). In order to avoid scratching polished surfaces and insulation, do not use abrasive cleaners.  
 Dry instrument thoroughly !

**Sterilisation**  
 May be autoclaved up to 137°C (280°F). Instruments must be completely exposed to the steam including inner surfaces.  
 Please follow the autoclave manufacturer's recommendations for handling, loading and parameter settings.  
 Do not sterilise with hot air.  
 Do not use 'Flash' autoclave procedures.

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