

5175 South Royal Atlanta Dr

Tucker, GA 30084 USA

770-414-4880

www.blackandblacksurgical.com

www.vorotekusa.com

INSTRUCTIONS FOR USE

DISPOSABLE MICROSUCTION TUBES AND CONTROLLER

WARNINGS:

- Black & Black Surgical instruments are to be used in accordance with these instructions for use. Read all sections of this insert prior to use. Improper use of the instrument may cause a serious injury.
- Black & Black Surgical instruments should be handled and operated by healthcare professionals completely familiar with their use, assembly, and disassembly. Use instruments for their intended purposes only.
- The Black & Black Surgical disposable microsuction tubes and controllers are packaged sterile and are intended for single use only. Re-use could result in infection or injury to the patient.

CAUTION: FDA (USA) law restricts this device to sale by or on the order of a physician.













INDICATIONS FOR USE:

These disposable microsuction tubes and controllers are indicated for cerumen management.

PACKAGED STERILE

INTENDED FOR SINGLE USE

REF	Gauge	Length	Hub Color
B78629	12 (2.7MM)	70MM	Blue
B78632	14 (2.0MM)	70MM	Green
B78637	17 (1.4MM)	70MM	Purple
B78624	Suction Controller		Black

	Reference number		Do not re-use
	Sterilized by Ethylene Oxide		Do not use if package is damaged
	Lot number		Keep dry
	Use by		Keep away from sunlight
	Manufacturing date		Do not re-sterilize
	Consult the instructions for use		Temperature limits