

ContourFoam™ Post-operative Liposuction Foam

Item # B89600 (Roll) and # B89601 (Cut Pieces)

Device Description:

ContourFoam™ is a flexible open-cell elastomeric polyurethane post operative foam designed to provide even, smooth padding and compression for added comfort and contouring. Our foam comes in a 6' x 14" x ½" roll. It can be cut to your desired shape and size, and then easily inserted between sterile dressing and compression garment. Packaged with UV-protected bags for longer shelf life.

Indications for Use:

ContourFoam™ is designed to be an external wrap to be used after a liposuction procedure. The main objective of the foam is to help reduce swelling and bruising after a liposuction procedure. It should be used according to doctors post operative instructions. On the IFU, there are no clinical applications set out or explained.

Warnings and Precautions:

Do not remove foam from bag until ready to use. Exposure to light for an extended period of time may cause foam to yellow. This yellowing does not affect the usefulness of the foam.

Do not use in the presence of flammable liquids! The foam is flammable when exposed to a flame at temperatures above 700-800 degrees F.

Failure to observe these cautions and contraindications may result in injury, malfunction or other unanticipated occurrences or events for the operator, staff and/or the patient.

Non-Sterile; Single Use Product