

## Nasal Tip Cartilage Punch System



*Designed to cut harvested cartilage into desired shapes, the Nasal Tip Cartilage Punch System features a manual reusable punch used for shaping cartilage during nasal surgery. The system saves time during preparation of nasal cartilage for re-introduction of nasal tip graft and provides more consistent results, as compared to free cutting operations using standard surgical instruments such as knives or scissors.*



*Designed by Bahman Guyuron, M.D.*

- B 67240 Punch, heart shape with handle
- B 67241 Punch, figure-eight shape with handle
- B 60760 Guyuron-Debaquey Ptosis Clamp
- B 67222 Guyuron Lateral Cruz Stabilizer Clamp

© 2019 Black & Black Surgical 25-0538 Rev A