

LipoPort™



LipoPort™ is a must-have for any surgeon performing liposuction. It improves the aesthetic outcome of incisions and aids in wound healing by saving the skin from friction burns. It can be used with both manual and power-assisted liposuction procedures. This product is reusable and made in the USA.

- SAVES THE SKIN FROM FRICTION BURNS AND AIDS IN WOUND HEALING
- IMPROVES AESTHETIC OUTCOME OF INCISIONS
- REUSABLE / AUTOCLAVABLE
- EASY TO INSERT AND ANCHOR WITH 3-0 NYLON SUTURE
- ACCOMMODATES STRAIGHT AND CURVED POWER-ASSISTED, LASER-ASSISTED, AND STANDARD CANNULAS
- EXTERNAL NOTCH FACILITATES OPTIMAL CANNULA ORIENTATION
- TWO CONVENIENT SIZES THAT ACCOMMODATE MOST CANNULAS

B 89LP34 Lipoport Medium accommodates 2-4mm cannula

B 89LP56 Lipoport Large accommodates 3-5mm basket cannula and up to 6mm standard cannula

LipoPort designed by John Connors, MD, plastic surgeon, Atlanta, GA.

© 2020 Black & Black Surgical | 25-0607 2019-8

Tucker, GA 30084
877.252.2517
770.414.4880
770.414.4879 Fax
www.blackandblacksurgical.com

b BLACK & BLACK
S U R G I C A L