

# MARKMANN DRAIN INSERTER



Frustrated with ineffective methods for placing drains after liposuction, Dr. Markmann developed this unique two-part drain inserter featuring a precision cannula and inserter — allowing the placement of the drain from the mons pubis liposuction site to the lower abdomen, as well as the waist into the flanks.



## Part Number B89610

- PRECISE DRAIN PLACEMENT REDUCES CHANCE OF SEROMAS •
- ALLOWS DRAIN PLACEMENT EVEN WITH CIRCUMFERENTIAL LIPOSUCTION •
- 36CM LENGTH AND 6.5MM ID ALLOWS FOR UP TO 15 FRENCH BLAKE™ OR SAFETY DRAINS •
- AUTOCLAVABLE, REUSABLE SYSTEM FOR MULTIPLE USES •

© 2018 Black & Black Surgical | 25-0607 2019-8

Tucker, GA 30084  
877.252.2517  
770.414.4880  
770.414.4879 Fax  
[www.blackandblacksurgical.com](http://www.blackandblacksurgical.com)

 **BLACK & BLACK**  
SURGICAL