

Infiltration Pump B89055

ULTRA FAST PERFORMANCE

FLUID COUNTING
FEATURE



VITRUVIAN™ Infiltration Pump B89055

The Black & Black Vitruvian™ Counting Infiltration Pump is a peristaltic pump that uses three rollers to create a smooth pumping action to deliver tumescent fluid to the operative site. Powered by a state of the art stepper motor, this machine provides a smooth and consistent pumping action at a rate of up to 999mL/min, making it the fastest counting pump in the market.

The Black & Black Vitruvian™ Infiltration Pump was designed with industry leading features such as fluid counting technology, a rate adjusting knob, prime button, reverse switch, and dual foot pedals that can be used in momentary or maintained mode. Large, bright LEDs are easily viewed from anywhere in the room and clearly display both volume and rate.

For best performance use Black & Black Softouch Pump Tubing (B89213).



- Peristaltic Pump with three rollers
- Momentary/ maintained pneumatic foot-switch functions for ON/OFF action or continuous action, respectively
- Dual voltage motors 100V/240V
- Bright, easy-to-read display screen
- Volume (mL) and Rate (mL/min) displayed
- Reverse switch for flow direction preference
- Prime button for fast filling of tubing line
- Dual foot pedals

B89055	VITRUVIAN™ INFILTRATION PUMP
B89213	13' SOFTOUCH PUMP TUBING SET-SINGLE SPIKE, BOX OF 10



5175 South Royal Atlanta Drive
 Tucker, GA 30084
 770.414.4880 phone
 1.877.252.2517 toll free
 770.414.4879 fax
 info@blackandblacksurgical.com
 www.blackandblacksurgical.com